

The New You

ULTRAFORMER III is the world's first and only Micro and Macro Focused Ultrasound (MMFU) device for high-intensity skin lifting and tightening. It is the world's most advanced and comprehensive system with dual handpieces and 7 interchangeable multi-depth cartridges.

EXPORT ONLY
EXPORT ONLY
UF3B012022



SYSTEM SPECIFICATIONS

Energy Type	HIFU (High Intensity Focused Ultrasound)
Power	0.1 ~ 3.0J (0.1J/step)
Spacing	1 ~ 2mm (0.1mm/step)
Length	0.5 ~ 25mm (5mm/step)
Rating	100-240 V~, 50 / 60Hz
Weight	35Kg
Dimension	500 x 515 x 1310mm
Cartridge	L7-1.5 (7MHZ/1.5mm) L7-3.0 (7MHZ/3.0mm) L4-4.5 (4MHZ/4.5mm) MF2 (5.5MHZ/2.0mm) MF6 (2MHZ/6.0mm) MF9 (2MHZ/9.0mm) MF13 (2MHZ/13.0mm)



T. +82-2-517-2114 | classys.com | ultraformer.com | info@classys.com
Copyright to CLASSYS®2022
Images and texts are intellectual property of Classsys.
Copying of this material can be subject to charges of both civil and criminal law of legal justice.



THE
NEW
YOU

Face Lifting +
Skin Tightening +
Face & Body Contouring +

ULTRAFORMER III

MMFU | Micro & Macro Focused Ultrasound

FACE LIFTING, SKIN TIGHTENING, and FACE & BODY CONTOURING

ULTRAFORMER III is an optimum, non-invasive ultrasound device that lifts and tightens the face for a youthful complexion and tightens the body for slimmer contours.

With precision per shot pulse, HIFU-powered transducers are designed to either remodel collagen to remove facial wrinkles and sagging skin or tighten body tissues to reflect your full potential.



Why Doctors Say That ULTRAFORMER III is the Answer



Micro & Macro Focused Ultrasound Technology

We comply with both clinical and international regulatory standards to assure safe and effective procedures with MMFU technology.



Multi-Function & Multi-Depth Cartridges

Delivery of an equal and fast distribution of ultrasound waves into multiple layers of skin and body tissues without direct contact to the epidermis.



Non-Invasive Procedure & Less Pain

Ultrasound energy equipped multi-depth cartridges of the ULTRAFORMER III only target specific layers of skin without damaging nearby tissue.



Fast & Precise Shot Application

HIFU-powered transducers are designed to either remodel collagen to remove facial wrinkles and sagging skin or tighten body tissues.



Dual Engine : High Peak Power Transducers

High peak power delivers energy effectively to target layers of skin without diffusing heat to surrounding tissues.



7 Interchangeable Cartridges for Face & Body

Unique multi-functional cartridges assisted by MMFU Technology are each designed for specific face and body applications.



Easy, User-Friendly Design and High Usability

ULTRAFORMER III has an interactive interface where the user can easily understand and navigate through the application in an efficient way.


HIGHLY ADVANCED MICRO & MACRO ULTRASOUND TECHNOLOGY MECHANISM OF MMFU

HIFU transducers deliver an equal distribution of ultrasound waves into multiple layers of skin and body tissues without direct contact with the epidermis.

Fitted with multiple cartridges that generate thermal energy, the ULTRAFORMER III maintains stabilized temperatures of 65–75 °C to induce coagulation for collagen renewal in the face and tightens body tissues for contouring waistlines, the thighs, and other parts of the body.


MICRO FOCUSED

With highly concentrated ultrasound beams, cartridges deliver power with accuracy and precision to the specified focal region.




MACRO FOCUSED

Cartridges with lower frequencies allow deeper penetration to the focal region with a larger energy density.




Dual Engine High Peak Power

High peak power delivers energy effectively to target layers of skin without diffusing heat to surrounding tissues, yielding less pain and a low rate of side effects combined with effective clinical results.




MMFU Thermal Energy Delivery

Generating approximately 65.4 °C into target layers induces natural coagulation and wound-healing processes to stimulate collagen renewal underneath the epidermis.



With 7 Interchangeable Cartridges FAST & PRECISE SHOT APPLICATION



PROVEN EFFICACY THROUGH CLINICAL STUDIES

Download from www.ultraformer.com/marketing-kit/

TIGHTENING EFFECTS OF HIGH-INTENSITY FOCUSED ULTRASOUND On Body Skin And Subdermal Tissue

Origin : [J EurAcadDermatolVenereol. 2016 Sep;30(9):1599-602.]

Procedures with the ULTRAFORMER III to increase skin elasticity in the face and tighten areas of the body through the non-invasive method of HIFU under multiple depth cartridges were evaluated to determine the efficacy of the device as a practical alternative to invasive procedures and surgery. The study was based on the Global Aesthetic Improvement Scale (GAIS) standards, and elasticity changes in the skin were measured via a cutometer.


For body tightening, we applied HIFU using transducers with a lower frequency and deep focal depth to effectively deliver ultrasound energy to skin tissues. HIFU appears to be a safe and effective treatment modality for dermal and subdermal tightening.

EFFICACY & SAFETY OF NON-INVASIVE BODY TIGHTENING With High-Intensity Focused Ultrasound (HIFU)

Origin : [Skin Research and Technology 2017; 1-5]

Clinical researchers have conducted treatment procedures on subjects with the ULTRAFORMER III to determine whether high-intensity focused ultrasound (HIFU) is an effective, non-invasive alternative to incision-based procedures or surgeries for improving skin elasticity and contours of the face and body. Over 12 weeks following treatment, researchers evaluated the device for its efficacy and safety based on the standards of the Global Aesthetic Improvement Scale (GAIS) and changes to skin elasticity measured over time via a cutometer.

* SG AIS* Improvement Rate _ Graphs for 12 Weeks Post Procedure



In this study, we observed significant improvements in two body regions (abdomen and thighs) as well as the cheek when targeted for HIFU treatment. By using newly developed transducers with different energy outputs and focal depths, HIFU treatment can be tailored to meet the unique physical characteristics of each patient.

* SG AIS : Subject Global Aesthetic Improvement Scale

* IG AIS : Investigator Global Aesthetic Improvement Scale

Clinical Studies from All Over the World

High-Intensity Focused Ultrasound: A Satisfactory, Non-invasive Procedure for Crow's Feet Wrinkles
Gyu Sik Jung et al. / South Korea

Efficacy of High Intensity Focused Ultrasound (HIFU) for Lifting and Tightening Lax Facial & Neck Skin
Sharmila Nayak / India

A New Treatment Protocol of Micro-Focused Ultrasound for Lower Eyelid Fat Bulging
Jeon Hye-Chan et al. / South Korea

Efficacy and Safety of Non-invasive Body Tightening with High Intensity Focused Ultrasound (HIFU)
Ko E.J. et al. / South Korea

The efficacy of macro-focused ultrasound in the treatment of upper facial laxity: A pilot study
Rungsima Wanitphakdeedecha, MD / Thailand

A study of efficacy and safety of high-intensity focused ultrasound for the treatment of melasma in Asians:
Vasanop Vachiramon, MD / Thailand

Non-Invasive Arm Fat Reduction
Logan William Thomas et al. / USA

High Speed Low-pain Micro Focused Ultrasound Tightening of the Lower Face and Neck
Adrian Lim, MD / Australia

SMAS Face Lift with HIFU Technology (High Intensity Focused Ultrasound) for the ULTRAFORMER Unit
Klaus Fritz, MD / Germany

Evaluation of Micro Focused Ultrasound for Lifting and Tightening the Face
Lee In-Ho et al. / South Korea

The Most Exciting International Evolution in the Non-surgical Facelift
Serena Lim, MD / Australia

Face Lifting and Body Modeling without a Scalpel
Radoslaw Rzepnikowski, MD / Poland

Tightening Effects of High Intensity Focused Ultrasound on Body Skin and Subdermal Tissue: A Pilot Study
Choi S.Y. et al. / South Korea



CLINICAL RESULTS AFTER

FACE LIFTING, SKIN TIGHTENING & BODY CONTOURING

FACE | Lifting & Tightening



Ref. Individual results may vary. All are unretouched photographs.

With ULTRAFORMER III

After Fast & Precise Shot Application

Participants shown in the following clinical photos received one ULTRAFORMER III treatment on the target areas presented. Evidence of clinical photos obtained from baseline and post-procedure measurements indicates results in progress.

Tightening & Contouring | BODY



TESTIMONIALS

by Medical Professionals



Dr Benjamin Ascher / MD. Plastic Surgeon

"I was deeply impressed by the skin-lifting results."

"I was deeply impressed by the skin-lifting results after experiencing ULTRAFORMER III for the first time. It shows significant scientific improvements of HIFU technology to develop great results. I believe CLASSYS will be one of the leading global companies to come based on their outstanding technology."



Dr Jean Louis Sebagh / MD. Plastic Surgeon

"ULTRAFORMER III is a game changer."

"I've been using the ULTRAFORMER III for more than a year now and I am very satisfied with its clinical results. I believe the ULTRAFORMER III is a game changer for skin lifting. Injections and surgery were our only options for skin lifting prior to HIFU, but now, the ULTRAFORMER III reflects better as a non-surgical treatment."



Dr Nobuhiro Suetake / MD. Plastic Surgeon

"So precise and brings immediate and gradual results."

"With multiple cartridges, you are able to customize treatments based on the patient's body shape and fat layer. It is so precise and brings immediate and gradual results. It's very impressive."



Dr Klaus Fritz / MD. Plastic Surgeon

"Painless treatment with high efficacy and no side effects."

"ULTRAFORMER III brings immediate tightening effects with different cartridges for the face. Surgery was the normal option for skin lifting procedures but now we can achieve equal results using ULTRAFORMER III non-invasively. My patients are always surprised about their remarkable clinical outcomes."



Dr Gerd G. Cauglitz / MD. Dermatologist

"Impressed by the results and have satisfied patients."

"We have had the ULTRAFORMER III now for quite a while. We are impressed by the results and have satisfied patients. From my point of view, it is an impressive and effective treatment, approaching a phase in a multilayer way, and involves having no downtime with results surfacing surprisingly quickly after the first procedure."



Dr Kim Jong-Seo / MD. Plastic Surgeon

"More effective and less painful."

"I am very impressed with ULTRAFORMER III in which there is a high rate of patient satisfaction. HIFU devices have replaced the trend of RF devices recently because these devices, like ULTRAFORMER, are more effective and less painful. It helps with blood circulation and remodels collagen effectively."



Dr Katia Paskova / MD. Dermatologist

"Suitable for patients who do not want surgical facelift procedures."

"The effect of ULTRAFORMER III is instantaneous. This is because we work on the tissues and various layers of the skin. With age, skincare products and cosmetics are just not enough. ULTRAFORMER is suitable for patients who do not want surgical facelift procedures."



Dr Marcin Ambroziak / MD. Dermatologist

"Very versatile device for sagging skin of the face and body."

"I have been using ULTRAFORMER III for more than 6 months now. This is a very versatile device because there are cartridges ranging from 1.5mm up to 9.0mm. So, we can use the device not only to treat sagging skin on the face, but we are able to treat other areas of the body as well."



Dr Michael H. Gold / MD. Dermatologist

"Patients are so satisfied with the great results reflected."

"For those who have never tried the ULTRAFORMER III with HIFU technology, I'd like to say that this is an outstanding tightening procedure device. We prefer to use it for tightening from the eyes down to the neck. Patients are so satisfied with their great results."



Dr Pedro Alberto Santos / MD. Plastic Surgeon

"Highest rate of satisfaction to both doctors & patients."

"I'm sure that the ULTRAFORMER III is the best device for providing the highest rate of satisfaction to both doctors and patients among other HIFU devices in the world. It's lucky for me to have a chance to use ULTRAFORMER III."



Dr Ishikawa Hirokazu / MD.

"Superior treatment capabilities over any other device."

"I've utilized many of the HIFU products available from all over the world. However, I can say that the ULTRAFORMER III is the best given its superior and more satisfying treatment capabilities over any other device."



Dr Franco Lauro / MD. Plastic Surgeon

"Skin tightening effects are amazing."

"The ULTRAFORMER III's skin tightening effects are amazing. Laser technology has limitations in terms of clinical results and side effects. The ULTRAFORMER III presents excellent skin tightening using HIFU that delivers energy, especially to the SMAS layer. The HIFU energy by ULTRAFORMER III takes apart old collagen and induces the formation of new collagen effectively."



Dr Atchima Suwanchinda / MD.

"Most satisfying treatment results for a couple of years."

"I'm sure that, no doubt, the ULTRAFORMER III is the best device for providing the happiest clients and most satisfying treatment results for a couple of years, even without any break-down. It is a great experience for me to have such an outstanding result."



Dr Rungsima Wanitphakdeedecha / MD. Plastic Surgeon

"So happily satisfied with immediate results."

"I serve customized treatments for both the face and the body with the various ULTRAFORMER III cartridges. Customers are so happily satisfied with immediate results. The ULTRAFORMER III is the most satisfying and best device for me."



Dr Adrian Lim / MD. Dermatologist

"Safe, effective, high-speed, low-pain procedure."

"MMFU therapy with the ULTRAFORMER III is a safe, effective high-speed, low-pain procedure that meets a clear need amongst patients seeking skin tightening."



Dr Felix Kusmana / MD. Biomed (AAM)

"Achieving the best results all over the body."

"We all know ULTRAFORMER III is the best HIFU machine that we must have in the clinic. Patients can achieve the best result almost all over their body."



Dr Jonathan Te-Peng Tseng / MD. Dermatologist

"So satisfying with a variety of treatments and services."

"I'm serving a variety of treatments using ULTRAFORMER III. It is so satisfying and interesting for me."



Dr Nedhal Khalifa / MD.

"Match patient tolerance and treat small areas preciously."

"It's an amazing fact for me that ULTRAFORMER III provides unique control that allows the user to match patient tolerance and treat small areas preciously."



Dr Dany Touma / MD.

"Quickly noticeable results and tolerable sessions."

"Definitely, ULTRAFORMER III differentiates itself from other HIFUs by its quickly noticeable results, its tolerable sessions, and having no adverse reactions."



Dr Irina Rakhova / MD.

"Makes my patients and me very very happy."

"I gave 2 nicknames to ULTRAFORMER III : Michelangelo's knife and a Time Machine. ULTRAFORMER III brings back youth rather than creating a new face with surgery. Each layer of skin regenerates collagen and gives your face a fantastic rejuvenation without any surgery or recovery time. ULTRAFORMER III makes my patients and me very very happy."